

Electronic dialysis Chair

A-GB-DIAC-70

1 Humanization design

► Modern design, high quality materials

Modern detailed finishing of the upholstery as well as mechanical and electrical parts underline the high standard of the comfort-Line. Continuous adjustment to all sitting and lying positions, adjustable armrest, utmost comfort for your patients. It makes the chairs suitable for dialysis centers, blood bank, blood donor centers and all kind of luxury hospital rooms.




► Multi-position adjustment

Option of choosing three **Denmark** motors for maximum patient comfort and reduction in staff work load. Backrest, seatrest, legrest, height are fully adjustable.

2 Saving space

Only 1.2Sq. Mtr. space required, saving about 40% space than the tradition medical bed.



3 Dust-proof & easily disinfect design

Good ABS covering for protecting the chair not corrosion, waterproof, dustproof, easily clean & disinfect & maintain.

4 High quality products

Taking the aero seat as the prototype improvement, upholstery with polyurethane, luxurious, comfortable and durability. And warm & light color makes patient enjoy the treatment.



▲ hands control



▲ Adjustable armrest



▲ adjustable footrest



▲ ABS covering



▲ castor with brake

Optional items

Tray or table	800X400X750/1050mm
I.V.pole	Standard I.V rope
TV	China famous brand 15"
Screen	Two or three fold

Technical data

Seat width	600mm	Back adjust	-12° - 75°
Total length	1900mm	Leg adjust	-80° - 12°
Seat width incl.arms	880mm	Trendelenburg	-12°
Height	560-820mm	Covering	ABS
Caster	100mm	Motors	3
Cushion	PU leather	Power supply	Ac110-240V 50/60Hz
Safe maximum load	240kg	Total weight	90Kg

Optional Color: green ■ blue ■ brown ■