

HOSPITAL BED SERIES

A-GB-care-001



ACOTA

Verweylaan 10,9725GM,
Groningen, The Netherlands
Email: info@acota-medical.com
www.acota-medical.com

ACOTA B.V.

ACOTA



Features

With the latest caring experiences and technology from Europe and USA, Acota B.V. From Netherlands brings you the latest model A-GB-care series.

It is an intelligent delivery bed designed based on the human engineering. It can be used for gynecology, delivery and recovery stage or be used for gynecological examining, diagnosing and operation.

The bed pays more attention t the emotion of the mother and baby. The convenient and simple operation method ultimately reduces the working stress of the doctors and nurses.

The operation is simple, flexible and convenient. The central brake systems makes the bed more stable.

Clear marked control panels are setting on the guardrail on the both sides where both nurses and patients can easily reach to.

Blow molding plastic side rails give the warm and comfortable feelings to the nursing staff and the patients.

Medical blowing plastic technology, hidden assistant table made from mould secures the good quality to the users

The retractable guiderail allows the bed to get through the narrow aisle more easily and realize the Zero-spare transportation to the crash truck.

The anti-bumpers where located on the head section provide good transportation method during the delivery from hospitals.

Guardrail is equipped with inductive light which much help the nursing staff and patients during the night and emergency transportation.

The UPS back-up batteries offer emergency use even the main power is off

The low bed position gives the patients more stable support.

With various labor positions (e.g. supine position, sitting position, lateral position, kneeling-squatting position), comforts different lying-in women.

The leg holder is retractable and height adjustable, which much more reduces the stress of the nursing staff.

The auxiliary table board offers a larger working space.

With the easy-operation foot rest and CPR spanner gives the best controlling feelings to the uses.

With the medical molding plastic mattress support, the features of anti-corrosion, anti-door , easy-clean can be realized easily.

Memory foam mattress gives the best comfort to the patients during all the positions.



Technical data

Bed dimension including bed and bumper	2350x940mm LxW
Net bed dimension	1950x880mm LxW
Bed height	Lowest ≤550mm , highest ≥890mm (exclude mattress)
Back section turning	≥65 degree
Seat section turning	≥12degree
Trendelenburg	≥8 degree
Leg section lifting distance	170mm
Footrest swing out	≥ 90degree
Footrest folding	≥ 90degree
Guardrail up/down movement	405mm

Standard delivery

A-GB-care-001	1 unit
IV stand	1 piece
Filth basin	1 piece
Mattress	1 piece