

ACOTA

LOOK NO POWER LEADS IS AVAILABLE!

AC Power & C-Cells batteries Power Diagnostic Unit

The unit is combined with two ways power supply, AC power and Batteries power supply, with the batteries power supply, it is Wireless unit where the user can place in any place.



Diagnostic Unit



The DIA wireless series provide all the benefits of traditional conventional wall units such as instruments security and light output, but also can be used as wireless use such as combination with a trolley.

- Brilliant illumination from the TrueTone LED instruments
- LED modules guaranteed for over 10 years!
- Over 6 months use from 3x C cells (depends on frequency and time to use)
- Interchangeable LED bulbs brings brilliant true illumination
- Simultaneous handle operation with no reduction in illumination
- Superb build quality, easy to clean and no sharp edges
- Secure handle cradles. Designed to be best available, no more accidental dislodging of handles from cradles with possible instrument head damage
- Auto On/Auto Off, Even if cross docked
- Out of cradle Auto Off after 5 minutes
- Easy to replace batteries, uses conventional alkaline C Cells, easily available and inexpensive
- Ideal for trolley mounting with no mains cable trip hazard






Our Version is for total quality

LOOK
NO POWER LEADS IS AVAILABLE!

DIA 50 Diagnostic Unit

LED Coaxial Ophthalmoscope
LED Fiber Otoscope
2 x 12ft spiral cord handles with Auto On / Auto Off
Cradle
Dispenser



Coaxial Ophthalmoscope	Fiber Otoscope	Cradle C-cell battery	Cradle C-cell battery & AC power	Dispenser
				
0010	0010.0	0100.1	0100.2	0060

Please specify whether Welch Allyn or Heine Otoscope tip fitment is required. Welch Allyn tip fitment supplied as standard, Heine upon request at no extra cost. Dispenser standard includes 2.7mm 35 pieces, 4.2mm 35 pieces.

DIA 50

Diagnostic set includes: 0010, 0010.0, 0100.2, 0600

DIA 50.1

Diagnostic set includes: 0010, 0010.0, 0100.1, 0600

DIA.50.T Mobile Trolley unit

ACOTA B.V.

Verweylaan 10,9725GM,
Groningen, The Netherlands
www.acota-medical.com

